

Anti-KAT6B antibody (400-500) (STJ119314)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 231kDa Observed Mw: 280kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23522
Gene Symbol	KAT6B
UniProt ID	KAT6B_HUMAN
Immunogen Region	400-500
Immunogen Sequence	DLDFVKQAQELSWKEICES GVEDCGRYPSVIEFGKYEIQ TWYSSPYPQEYARLPKLYLC EFCLKYMKSKNILLRHSSKKG GWFHPPANEIYRRKDLSVFE V
Specificity	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human KAT6B (NP_001243398.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using KAT6B Rabbit polyclonal antibody (STJ119314) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.