


Anti-ASB4 antibody (1-300) (STJ119327)

GENERAL INFORMATION		
Product Type	Primary antibodies	
Applications	IHC-P/ELISA	
Host / Source	Rabbit	
Reactivity	Human/Mouse	


PRODUCT PROPERTIES		
Clonality	Polyclonal	
Concentration	Lot specific	
Conjugation	Unconjugated	
Purification	Affinity purification	
Dilution Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.	
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.	
Isotype	IgG	
Molecular Weight	Protein Mw: 48kDa Observed Mw: Refer to figures	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.	

TARGET INFORMATION		
Gene ID	51666	
Gene Symbol	ASB4	
UniProt ID	ASB4_HUMAN	
Immunogen Region	1-300	
Immunogen Sequence	MDGTTAPVTKSGAAKLVKRN FLEALKSNDFGKLKAILQR QIDVDTVFEVEDENMVLASY KQGYWLPSPYKLSWATGLH LSVLFQGHVECLLVLDHNAT INCRPNKGKTPLVHACEMANV DCVKILCDRGAKLNCYSLSG HTALHFCTTPSSILCAKQLV WRGANVNMKTNNQDEETPLH TAAHFGLSELVAFYVEHGAI VDSVNAHMETPLAIAAYWAL RFKEQEYSTEHLVCRMLL	
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ASB4 (NP_665879.1).	


ADDITIONAL INFORMATION		
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.	



Immunohistochemistry analysis of ASB4 in paraffin-embedded mouse kidney using ASB4 Rabbit polyclonal antibody (STJ119327) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of ASB4 in paraffin-embedded human uterine cancer using ASB4 Rabbit polyclonal antibody (STJ119327) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of ASB4 in paraffin-embedded human liver using ASB4 Rabbit polyclonal antibody (STJ119327) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.