

Anti-NCOA3 antibody (1-100) (STJ119458)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 155kDa Observed Mw: 160kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	8202
Gene Symbol	NCOA3
UniProt ID	NCOA3_HUMAN
Immunogen Region	1-100
Immunogen Sequence	MSG LGENLDPLASDSRKRL PCDTPGQGLTCSGEKRRREQ ESKYIEELAE LISANLSDID NFNVKPDKCAILKETVRQIR QIKEQGKTISNDDDVQKADV
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NCOA3 (NP_858045.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of HeLa, using NCOA3 rabbit polyclonal antibody (STJ119458) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunoprecipitation analysis of 300 Mu g extracts of HeLa cells using 3 Mu g NCOA3 rabbit polyclonal antibody (STJ119458). Western blot was performed from the immunoprecipitate using NCOA3 rabbit polyclonal antibody (STJ119458) at a dilution of 1:2000.