

Anti-VPS37B antibody (60-170) (STJ119806)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 31kDa Observed Mw: 31kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79720
Gene Symbol	VPS37B
UniProt ID	VPS37B_HUMAN
Immunogen Region	60-170
Immunogen Sequence	LLYQPQLDTLTKARLTQKYQE LQVLFEAYQIKKTKLDRQSS SASLETLLALLQAEGAKIEE DTENMAEKFLDGELPLDSFI DVYQSKRKLAHMRRVKIEKL QEMVLKGRQLP
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 60-170 of human VPS37B (NP_078943.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from Rat thymus, using VPS37B Rabbit polyclonal antibody (STJ119806) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.