

## Anti-PXDN antibody (1320-1479) (STJ119915)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                        |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                             |
| Concentration       | Lot specific                                                                                                                                           |
| Conjugation         | Unconjugated                                                                                                                                           |
| Purification        | Affinity purification                                                                                                                                  |
| Dilution Range      | WB:1:500-1:1000<br>ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation         | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                                       |
| Isotype             | IgG                                                                                                                                                    |
| Molecular Weight    | Protein Mw: 165kDa<br>Observed Mw: 165kDa                                                                                                              |
| Storage Instruction | Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                          |

### TARGET INFORMATION


|                    |                                                                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">7837</a>                                                                                                                                                           |
| Gene Symbol        | <a href="#">PXDN</a>                                                                                                                                                           |
| UniProt ID         | <a href="#">PXDN_HUMAN</a>                                                                                                                                                     |
| Immunogen Region   | 1320-1479                                                                                                                                                                      |
| Immunogen Sequence | RTRGQFNAFSYHFRGRRSL EFSYQEDKPTKTRPRKIPSV GRQGEHLSNSTSAFSTRSDA<br>SGTNDRFREFVLEMQKTITDL RTQIKKLESRLSTTECV DAG GESHANNTKWKKDACTICEC<br>KDGQVTCFVEACPPATCAVP VNIPGACCPVCLQKRAEEKP |
| Specificity        | Recombinant fusion protein containing a sequence corresponding to amino acids 1320-1479 of human PXDN (NP_036425.1).                                                           |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of Raji cells, using PXDN antibody (STJ119915) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Western blot analysis of extracts of Mouse lung, using PXDN antibody (STJ119915) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.