

Anti-Mature-BDNF antibody (C-Term) (STJ13100064)

GENERAL INFORMATION

| | |
|---------------|--------------------------|
| Product Type | Primary antibodies |
| Applications | IHC/WB |
| Host / Source | Guinea Pig |
| Reactivity | Human/Rat/Mouse/Marmoset |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Polyclonal |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | IHC, WB. A dilution of 1:500 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications. |
| Formulation | Lyophilised |
| Storage Instruction | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles. |

TARGET INFORMATION

| | |
|------------------|----------------------------|
| Gene ID | 627 |
| Gene Symbol | BDNF |
| UniProt ID | BDNF_HUMAN |
| Immunogen Region | C-Term |
| Specificity | Specific for mature BDNF. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|