

Anti-TRPC6 antibody (10-60 aa) (STJ13100423)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	WB. A concentration of 3-10 μ g/ml is recommended. The optimal concentration should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	22068
Gene Symbol	Trpc6
UniProt ID	TRPC6_MOUSE
Immunogen Region	10-60 aa
Specificity	Specific for TRPC6.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



WB on mouse brain lysate. Blocking: 1% LFDM for 30 min at RT; primary antibody dilution 5 μ g/ml incubated at 4°C overnight.