

Anti-OMP antibody (STJ13100485)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB
Host / Source	Guinea Pig
Reactivity	Human/Rat/Mouse/Marmoset/Hamster

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB. Use at a dilution of 1:500. The optimal dilution should be determined by the end user. For IHC, Zamboni or 4% PFA fixatives are recommended.
Formulation	Lyophilised
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	4975
Gene Symbol	OMP
UniProt ID	OMP_HUMAN
Specificity	Highly specific for mature olfactory neurons (including axon and terminals).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---