

## Anti-TGN38 antibody (STJ13100507)

### GENERAL INFORMATION

|               |                          |
|---------------|--------------------------|
| Product Type  | Primary antibodies       |
| Applications  | IHC/WB                   |
| Host / Source | Guinea Pig               |
| Reactivity    | Rat/Mouse/Human/Marmoset |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                                                                                                        |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                      |
| Purification        | Affinity purification                                                                                                                                                                                                                             |
| Dilution Range      | IHC, WB. A dilution of 1: 500 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.                                                                                                    |
| Formulation         | Lyophilised                                                                                                                                                                                                                                       |
| Storage Instruction | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles. |

### TARGET INFORMATION

|             |                                                       |
|-------------|-------------------------------------------------------|
| Gene ID     | <a href="#">22135</a>                                 |
| Gene Symbol | <a href="#">Tgoln2</a>                                |
| UniProt ID  | <a href="#">TGON2_MOUSE</a>                           |
| Specificity | Specific for TGN38 (or TGN46, TGN48, TGN51 in human). |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|