

**Anti-ABCA4 antibody (1450-1500 aa) (STJ13100654)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC/WB
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">24</a>
Gene Symbol	<a href="#">ABCA4</a>
UniProt ID	<a href="#">ABCA4_HUMAN</a>
Immunogen Region	1450-1500 aa
Specificity	Specific for ABCA4.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---