

Anti-Phospho-CRMP2-T514/S518 antibody (STJ13100794)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	IHC/WB
	Host / Source	Sheep
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. This antibody recognises mainly phospho T514 & S518 of CRMP-2 (but also it is quite likely to recognise the unphosphorylated form).
	Formulation	Lyophilised
	Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
TARGET INFORMATION		
	Specificity	Specific for CRMP2.
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.