

Anti-CRP antibody (Full Length) (STJ13101004)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Sheep
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA, IHC, WB. A dilution of 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	1401
Gene Symbol	CRP
UniProt ID	CRP_HUMAN
Immunogen Region	Full Length
Specificity	Specific for CRP.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---