

Anti-Olig2 antibody (STJ13101561)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB/FC
Host / Source	Guinea Pig
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB, Flow Cyt. A dilution of 1:500 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	50913
Gene Symbol	Olig2
UniProt ID	OLIG2_MOUSE
Specificity	Specific for Olig2.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---