

Anti-OR11H4 antibody (280-322 aa) (STJ13101720)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC. A concentration of 10-50 µg/ml is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	390442
Gene Symbol	OR11H4
UniProt ID	O11H4_HUMAN
Immunogen Region	280-322 aa
Specificity	Specific human olfactory receptor Olfactory receptor OR11H4.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---