

Anti-VN1A1, V1A13, V1A14 antibody (200-250 aa) (STJ13101850)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC/WB
Host / Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:2000 is recommended for WB and to 1:1000 for IHC-P. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Immunogen Region	200-250 aa
Specificity	Specific for V1A1, V1A13, V1A14.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---