

## Anti-TRPV5 antibody (4th Cytoplasmic Loop) (STJ13102367)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB
Host / Source	Sheep
Reactivity	Rabbit

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB. A dilution of 1:300 to 1:4000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

### TARGET INFORMATION

Immunogen Region	4th Cytoplasmic Loop
Specificity	Specific for TRPV5.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---