

Anti-EGFP/GFP/Venus antibody (STJ13102673)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB/ELISA/FC
Host / Source	Carp

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Dilution Range	Reconstitute in 500 ul of sterile water. Centrifuge to remove any insoluble material.
Formulation	Lyophilised. Contains 0.02% benzalkonium chloride.
Isotype	IgM
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20C for long term storage and refrigerated at 2-8C for a shorter term.

TARGET INFORMATION

Specificity	Specific for EGFP/GFP/Venus
-------------	-----------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---