

Anti-Rabbit Anti-mouse ImmunoSwitch antibody (STJ13102675)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | Immunoswitch |
| Host / Source | NZ white rabbit |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Clonality | Polyclonal |
| Conjugation | Unconjugated |
| Purification | Purified IgG |
| Dilution Range | Reconstitute in 500 ul DI water. |
| Formulation | Lyophilised. Contains 0.02% benzalkonium chloride. |
| Isotype | IgG |
| Storage Instruction | Maintain the lyophilised/reconstituted antibodies frozen at -20C for long term storage and refrigerated at 2-8C for a shorter term. |

TARGET INFORMATION

| | |
|-------------|------------------------|
| Specificity | Specific for mouse IgG |
|-------------|------------------------|

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|