

## Anti-Rabbit Anti-donkey ImmunoSwitch antibody (STJ13102680)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Immunoswitch
Host / Source	NZ white rabbit
Reactivity	Donkey

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Purified IgG
Dilution Range	Reconstitute in 500 ul DI water.
Formulation	Lyophilised. Contains 0.02% benzalkonium chloride.
Isotype	IgG
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20C for long term storage and refrigerated at 2-8C for a shorter term.

### TARGET INFORMATION

Specificity	Specific for donkey IgG
-------------	-------------------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---