

Anti-Guinea pig Anti-rabbit ImmunoSwitch antibody (STJ13102697)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | Immunoswitch |
| Host / Source | Guinea pig |
| Reactivity | Rabbit |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Polyclonal |
| Conjugation | Unconjugated |
| Purification | Purified IgG |
| Dilution Range | Reconstitute in 500 ul DI water. |
| Formulation | Lyophilised. Contains 0.02% benzalkonium chloride. |
| Isotype | IgG |
| Storage Instruction | Maintain the lyophilised/reconstituted antibodies frozen at -20C for long term storage and refrigerated at 2-8C for a shorter term. |

TARGET INFORMATION

| | |
|-------------|-------------------------|
| Specificity | Specific for rabbit IgG |
|-------------|-------------------------|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|