

Anti-Guinea pig Anti-pig ImmunoSwitch antibody (STJ13102700)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Immunoswitch
Host / Source	Guinea pig
Reactivity	Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Purified IgG
Dilution Range	Reconstitute in 500 ul DI water.
Formulation	Lyophilised. Contains 0.02% benzalkonium chloride.
Isotype	IgG
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20C for long term storage and refrigerated at 2-8C for a shorter term.

TARGET INFORMATION

Specificity	Specific for pig IgG
-------------	----------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---