

## Anti-Rat Anti-goat ImmunoSwitch antibody (STJ13102707)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Immunoswitch
Host / Source	Rat
Reactivity	Goat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Purified IgG
Dilution Range	Reconstitute in 500 ul DI water.
Formulation	Lyophilised. Contains 0.02% benzalkonium chloride.
Isotype	IgG
Storage Instruction	Maintain the lyophilised/reconstituted antibodies frozen at -20C for long term storage and refrigerated at 2-8C for a shorter term.

### TARGET INFORMATION

Specificity	Specific for goat IgG
-------------	-----------------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---