

## Anti-GFP antibody (STJ140005)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/IHC-P/IHC-F
Host / Source	Goat
Reactivity	GFP

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:5000 IF 1:50-1:1000 IHC-P 1:50-1:1000 IHC-F 1:50-1:1000 Cabaço LC Bento-Lopes L Neto MV et al. JID Innov 2022 Jun PMID: 36090299 Vyas P Wu JS Zimmerman A et al. J Assoc Res Otolaryngol 2016 Sep PMID: 27696081 Larson
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20 C for long-term storage. Store at 2-8 C for up to one month.

### TARGET INFORMATION


Immunogen Sequence	MVSKGEELFTGVVPILVELD VSKGEELFTGVVPILVELDG SKGEELFTGVVPILVELDGD KGEELFTGVVPILVELDGDV GEELFTGVVPILVELDGDVN EELFTGVVPILVELDGDVNG ELFTGVVPILVELDGDVNGH LFTGVVPILVELDGDVNGHK FTGVVPILVELDGDVNGHKF TGVVPILVELDGDVNGHKFS GVPILVELDGDVNGHKFSV VVPILVELDGDVNGHKFSV
Specificity	In 293HEK cells transfected with cds plasmid detects a band of 27 kDa by Western blot. This antibody does recognize YFP and CFP and does not recognize RFP and mCherry fluorescent proteins.

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Anti GFP Goat Polyclonal antibody immunofluorescence drosophila STJ140005



Anti GFP Goat Polyclonal antibody immunofluorescence STJ140005



Anti GFP Goat Polyclonal antibody immunofluorescence STJ140005



Anti GFP Goat Polyclonal antibody western blot hek293 STJ140005