

Anti-CLTC-H Chain antibody (50aa to N-Term) (STJ140022)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:250-1:1000 IF 1:25-1:250 Imbastari F PhD Thesis Humboldt-Universität zu Berlin Germany 2020 Cai Y Postnikova EN Bernbaum JG et al. J Virol. 2015 Jan 1:89 (1:844-56. PMID: 25355889
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Gene ID	1213
Gene Symbol	CLTC
UniProt ID	CLH1_HUMAN
Immunogen Region	50aa to N-Term
Immunogen Sequence	MAQILPIRFQEHLQLQNLGI AQILPIRFQEHLQLQNLGIN QILPIRFQEHLQLQNLGINP
Specificity	Detects endogenous levels of total CLTC by Western blot in the whole cell lysates (MDCK, At-T20, COS-7).

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti CLTC H Chain Goat Polyclonal antibody western blot STJ140022



Anti CLTC H Chain Goat Polyclonal antibody immunofluorescence STJ140022



Anti CLTC H Chain Goat Polyclonal antibody immunofluorescence STJ140022