

Anti-RPS6 antibody (190aa C-Term) (STJ140049)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000 IF 1:25-1:250
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Gene ID	6194
Gene Symbol	RPS6
UniProt ID	RS6_HUMAN
Immunogen Region	190aa C-Term
Immunogen Sequence	MALKKQRTKKNKEEAAEYAK ALKKQRTKKNKEEAAEYAKL LKKQRTKKNKEEAAEYAKLL
Specificity	Detects a band of approximately 30 kDa by Western blot in the following whole cell lysates: MDCK, AtT20, Jurkat and 3T3.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti RPS6 Goat Polyclonal antibody immunofluorescence STJ140049



Anti RPS6 Goat Polyclonal antibody western blot STJ140049