

## Anti-WEE1 antibody (1440aa to C-Term) (STJ140054)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Goat               |
| Reactivity    | U.maydis           |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                              |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                                   |
| Concentration       | 3 mg/mL                                                                                                                                                                      |
| Conjugation         | Unconjugated                                                                                                                                                                 |
| Purification        | This antibody is epitope-affinity purified from goat antiserum.                                                                                                              |
| Dilution Range      | WB 1:500-1:2000                                                                                                                                                              |
| Formulation         | PBS, 20% Glycerol and 0.05% Sodium Azide.                                                                                                                                    |
| Isotype             | IgG                                                                                                                                                                          |
| Storage Instruction | For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours. |

### TARGET INFORMATION

|                    |                                                                                                             |
|--------------------|-------------------------------------------------------------------------------------------------------------|
| Immunogen Region   | 1440aa to C-Term                                                                                            |
| Immunogen Sequence | MSARTSISTSMSSFSSHQAV SARTSISTSMSSFSSHQAVI ARTSISTSMSSFSSHQAVIR<br>RTSISTSMSSFSSHQAVIRV TSISTSMSSFSSHQAVIRVR |
| Specificity        | Detects the immunogen by Western blot.                                                                      |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|