

Anti-RAB1A antibody (110aa C-Term) (STJ140056)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/IHC-P/IHC-F
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:250-1:1000 IHC-F 1:250-1:1500 IHC-P 1:250-1:1500 IF 1:50-1:250 ELISA ND Mojallal M Zheng Y Hultin S et al. Nat Commun 2014 Aug 1:54557. PMID: 25080966
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Gene ID	5861
Gene Symbol	RAB1A
UniProt ID	RAB1A_HUMAN
Immunogen Region	110aa C-Term
Specificity	Detects total Rab1 protein by Western blot in the following human, rat and mouse whole cell lysates and transfected cells with GFP-Rab1a and GFP-Rab1b cds.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Anti RAB1A Goat Polyclonal antibody immunohistochemistry human STJ140056



Anti RAB1A Goat Polyclonal antibody immunohistochemistry mouse STJ140056



Anti RAB1A Goat Polyclonal antibody immunofluorescence cos7 STJ140056



Anti RAB1A Goat Polyclonal antibody immunohistochemistry human STJ140056