

Anti-RAB7A antibody (C-Term) (STJ140063)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/IHC-P/IHC-F
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:250-1:5000 IF 1:50-1:250 IHC-F 1:100-1:500 IHC-P 1:100-1:500 Cardoso MHS PhD Thesis NOVA University of Lisbon Portugal 2018
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in-20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Gene ID	7879
Gene Symbol	RAB7A
UniProt ID	RAB7A_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	FLIQASPRDPENFPFVVLGN LIQASPRDPENFPFVVLGNK IQASPRDPENFPFVVLGNKI QASPRDPENFPFVVLGNKID ASPRDPENFPFVVLGNKIDL SPRDPENFPFVVLGNKIDLE
Specificity	Detects Rab7a protein in the human, rat and mouse whole cell lysates and transfected cells with GFP-Rab7a by Western blot. This Ab is specific for Rab7a.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Anti RAB7A Goat Polyclonal antibody western blot
hela STJ140063




Anti RAB7A Goat Polyclonal antibody
immunohistochemistry human STJ140063



Anti RAB7A Goat Polyclonal antibody
immunohistochemistry human STJ140063



Anti RAB7A Goat Polyclonal antibody western blot
STJ140063



Anti RAB7A Goat Polyclonal antibody
immunofluorescence STJ140063