

Anti-RAB20 antibody (130aa C-Term) (STJ140072)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:250-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20 C for long-term storage. Store at 2-8 C for up to one month.

TARGET INFORMATION

Gene ID	55647
Gene Symbol	RAB20
UniProt ID	RAB20_HUMAN
Immunogen Region	130aa C-Term
Immunogen Sequence	MEGPASGKVGSCVSTKVPKQ EGPASGKVGSCVSTKVPKQV GPASGKVGSCVSTKVPKQVQ PASGKVGSCVSTKVPKQVQP ASGKVGSCVSTKVPKQVQPE SGKVGSCVSTKVPKQVQPED
Specificity	Detects Rab20 in transfected cells with GFP-Rab20 by Western blot. Detects levels of endogenous Rab20 (testis, liver and kidney) by Western blot.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti RAB20 Goat Polyclonal antibody western blot
kidney STJ140072