

Anti-RAB27A antibody (120aa C-Term) (STJ140074)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/IHC-P/IHC-F
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000 IHC-F 1:50-1:400 IHC-P 1:50-1:400 IF 1:50-1:250 ELISA ND Ferreira JV Soares AR Ramalho J et al. Sci Adv 2022 Mar PMID: 35333565 Ferreira JV Rosa Soares A Ramalho JS et al. PLoS One 2019 Oct. PMID: 31613922
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Gene ID	5873
Gene Symbol	RAB27A
UniProt ID	RB27A_HUMAN
Immunogen Region	120aa C-Term
Specificity	Detects total Rab27 levels by Western blot in the following human, rat and mouse whole cell lysates and transfected cells with GFP-Rab27a and GFP-Rab27b cds.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Anti RAB27A Goat Polyclonal antibody immunofluorescence rpe STJ140074




Anti RAB27A Goat Polyclonal antibody western blot STJ140074




Anti RAB27A Goat Polyclonal antibody immunofluorescence STJ140074



Anti RAB27A Goat Polyclonal antibody immunofluorescence rpe STJ140074



Anti RAB27A Goat Polyclonal antibody immunohistochemistry human STJ140074



Anti RAB27A Goat Polyclonal antibody immunofluorescence STJ140074