

Anti-RAB28 antibody (120aa C-Term) (STJ140079)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:250-1:2000 IF 1:25-1:250
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in-20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Gene ID	9364
Gene Symbol	RAB28
UniProt ID	RAB28_HUMAN
Immunogen Region	120aa C-Term
Immunogen Sequence	MSETQPLVALVGNKIDLEHM SETQPLVALVGNKIDLEHMR ETQPLVALVGNKIDLEHMRT TQPLVALVGNKIDLEHMRTV QPLVALVGNKIDLEHMRTVK PLVALVGNKIDLEHMRTVKA
Specificity	Detects Rab28 by Western blot in transfected cells with GFP-Rab28.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti RAB28 Goat Polyclonal antibody western blot hek293 STJ140079



Anti RAB28 Goat Polyclonal antibody immunofluorescence STJ140079