

Anti-RAB32 antibody (120aa C-Term) (STJ140080)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000 IF 1:100-1:500
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Gene ID	10981
Gene Symbol	RAB32
UniProt ID	RAB32_HUMAN
Immunogen Region	120aa C-Term
Immunogen Sequence	KNLDLSKVHLPNGSPIPAVL NDLDSKVHLPNGSPIPAVLL DLDSKVHLPNGSPIPAVLLA LDSKVHLPNGSPIPAVLLAN DSKVHLPNGSPIPAVLLANK SKVHLPNGSPIPAVLLANKC
Specificity	Detects an endogenous band by Western blot in Melan-ink, XB2 and 3T3 cell line lysates (approximately 30 kDa) and in transfected cells with GFP-Rab32 cds.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti RAB32 Goat Polyclonal antibody western blot
STJ140080



Anti RAB32 Goat Polyclonal antibody
immunofluorescence rpe STJ140080