

Anti-CAT antibody (STJ140091)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	847
Gene Symbol	CAT
UniProt ID	CATA_HUMAN
Immunogen Sequence	MADSRDPASDQMQRHWKEQRA ADSRDPASDQMQRHWKEQRAA DSRDPASDQMQRHWKEQRAAQ SRDPASDQMQRHWKEQRAAQK RDPASDQMQRHWKEQRAAQKA DPASDQMQRHWKEQRAAQKAD PASDQMQRHWKEQRAAQKADV ASDQMQRHWKEQRAAQKADVL SDQMQRHWKEQRAAQKADVLT DQMQRHWKEQRAAQKADVLT QMQRHWKEQRAAQKADVLTG MQRHWKEQRAAQKADVLTG
Specificity	Detects a band of 60 kDa by Western blot.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti CAT Goat Polyclonal antibody western blot
bovine liver STJ140091