

Anti-RAB35 antibody (110aa C-Term) (STJ140106)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:250-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	11021
Gene Symbol	RAB35
UniProt ID	RAB35_HUMAN
Immunogen Region	110aa C-Term
Immunogen Sequence	MNQNCDDVCRILVGKNKDDP NQNCDDVCRILVGKNKDDPE QNCDDVCRILVGKNKDDPER NCDDVCRILVGKNKDDPERK CDDVCRILVGKNKDDPERKV
Specificity	Detects Rab35 in cell lysates and transfected cells by Western blot.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti RAB35 Goat Polyclonal antibody western blot
hek293 STJ140106