

Anti-INS antibody (STJ140108)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	3630
Gene Symbol	INS
UniProt ID	INS_HUMAN
Immunogen Sequence	EAEDLQVGQVELGGGPGAGS AEDLQVGQVELGGGPGAGSL EDLQVGQVELGGGPGAGSLQ DLQVGQVELGGGPGAGSLQP
Specificity	Detects C-peptide using MBP-C-peptide recombinant fusion protein and GFP-C-peptide expressed in transfected mammalian cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti INS Goat Polyclonal antibody western blot
STJ140108