

Anti-mgp antibody (26aa C-Term) (STJ140111)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Zebrafish

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:250-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Region	26aa C-Term
Immunogen Sequence	ESRESLEVFVNPYQANAFMR SRESLEVFVNPYQANAFMRN RESLEVFVNPYQANAFMRNT ESLEVFVNPYQANAFMRNTQ
Specificity	Detects endogenous levels of MGP protein using 5 days zebrafish embryos lysates.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---