

Anti-TP53 antibody (280aa C-Term) (STJ140114)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:5000 IF 1:25-1:250 Alves LS Marques ARA Padrao N et al. J Cell Sci 2021 Sep PMID: 34528688
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Gene ID	7157
Gene Symbol	TP53
UniProt ID	P53_HUMAN
Immunogen Region	280aa C-Term
Immunogen Sequence	DRRTEEEENLRKKGEPHHELP RRTEEEENLRKKGEPHHELPP RTEEENLRKKGEPHHELPPG TEEEENLRKKGEPHHELPPGS EEENLRKKGEPHHELPPGST EENLRKKGEPHHELPPGSTK
Specificity	Detects endogenous levels of P53 by Western blot and in transfected cells with GFP-P53.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti TP53 Goat Polyclonal antibody immunofluorescence STJ140114



Anti TP53 Goat Polyclonal antibody western blot hela STJ140114