

Anti-GST antibody (STJ140117)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	GST

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000 Seixas C Choi SY Polgar N et al. Mol Biol Cell 2015 Nov 18. PMID: 26582389 Catarino S Ribeiro-Rodrigues TM Sa Ferreira R et al. Cells 2020 Apr. PMID: 32272685
Formulation	PBS only
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Sequence	MSPILGYWKIKGLVQPTRL L SPILGYWKIKGLVQPTRL L L PILGYWKIKGLVQPTRL L L E ILGYWKIKGLVQPTRL L L EY LGYWKIKGLVQPTRL L L EY L GYWKIKGLVQPTRL L L EY L E YWKIKGLVQPTRL L L EY L E E W K I K G L V Q P T R L L E Y L E E K K I K G L V Q P T R L L E Y L E E K Y IKGLVQPTRL L L EY L E E K Y E K G L V Q P T R L L E Y L E E K Y E E G L V Q P T R L L E Y L E E K Y E E
Specificity	In lysates of E. coli cells transformed with pGEX plasmid detects a band of 26 kDa by Western blot.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Anti GST Goat Polyclonal antibody western blot
STJ140117