

Anti-VFP antibody (STJ140120)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/IHC-P/IHC-F/IEM
Host / Source	Goat
Reactivity	YFP/GFP/EGFP/Venus

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000 IHC-F 1:50-1:1000 IHC-P 1:50-1:1000 IF 1:50-1:1000 IEM 1:50-1:1000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Specificity	In 293HEK cells transfected with cds plasmid detects a band of 27 kDa by Western blot. This antibody also recognizes Venus/GFP and does not cross-react to mCherry/red fluorescent proteins.
-------------	--

ADDITIONAL INFORMATION


Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Anti VFP Goat Polyclonal antibody western blot
hek293 STJ140120



Anti VFP Goat Polyclonal antibody
immunoelectron microscopy STJ140120



Anti VFP Goat Polyclonal antibody
immunofluorescence hek293 STJ140120