

Anti-Azurin antibody (STJ140126)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Bacteria/Pseudomonas/Burkholderia

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:5000 Unver Y Gun BS Acar M et al. Prep Biochem Biotechnol 2020 Dec PMID: 33346686 Yildiz M Askin H. Prep Biochem Biotechnol 2019 Jun. PMID: 31156029 PMID: 31156029 Paydarnia N Khoshtinat Nikkhai S Fakhravar A et al. Mol Biol Re
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Sequence	MLRKLLVPVASPLILLGPSS LRKLLVPVASPLILLGPSSA RKLLVPVASPLILLGPSSAE KLLVPVASPLILLGPSSAEC LLVPVASPLILLGPSSAEC LVPVASPLILLGPSSAECV VPVASPLILLGPSSAECSDV PVASPLILLGPSSAECSDI
Specificity	Detects azurin by Western blot.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Anti Azurin Goat Polyclonal antibody western blot
STJ140126