

Anti-ATXN3 antibody (120aa to 250aa) (STJ140128)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:250-1:2000
Formulation	PBS only
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	4287
Gene Symbol	ATXN3
UniProt ID	ATX3_HUMAN
Immunogen Region	120aa to 250aa
Immunogen Sequence	KEHWFTVRKLGKQWFNLNSL EHWFTVRKLGKQWFNLNSLL HWFTVRKLGKQWFNLNSLLT WFTVRKLGKQWFNLNSLLTG FTVRKLGKQWFNLNSLLTGP TVRKLGKQWFNLNSLLTGPE VRKLGKQWFNLNSLLTGPEL
Specificity	Detects levels of ATXN3 protein in transfected cells with GFP-ATXN3 and endogenous levels in U118, SH-SY5Y, At-T20 and hTERT cell lysates by Western blot.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Anti ATXN3 Goat Polyclonal antibody western blot
STJ140128