

Anti-VSV G-Protein antibody (501-511) (STJ140131)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	VSV-G-Tag

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:5000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Region	501-511
Immunogen Sequence	MYTDIEMNRLGKNPMYTDIE YTDIEMNRLGKNPMYTDIEM
Specificity	This antibody recognizes recombinant proteins containing YTDIEMNRLGK epitope tag fused to either amino- or carboxy-terminal of targeted proteins in transfected cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti VSV G-Protein Goat Polyclonal antibody
western blot STJ140131