

Anti-RAB9A antibody (115aa C-Term) (STJ140136)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:250-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	9367
Gene Symbol	RAB9A
UniProt ID	RAB9A_HUMAN
Immunogen Region	115aa C-Term
Immunogen Sequence	QNLSNWKKEFIYYADVKEPE NLSNWKKEFIYYADVKEPES LSNWKKEFIYYADVKEPESF SNWKKEFIYYADVKEPESFP NWKKEFIYYADVKEPESFPF WKKEFIYYADVKEPESFPFV
Specificity	Rab9a antibody detects endogenous levels of Rab9a protein and transfected cells with GFP-Rab9a cds.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti RAB9A Goat Polyclonal antibody western blot
STJ140136