

Anti-RAB1B antibody (110aa C-Term) (STJ140139)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IHC-F
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Bovine/Canine/Chicken/Avian/Donkey/Feline/Goat/Guinea Pig/Hamster/Horse/Porcine/Rabbit/Sheep/Simian

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:5000 IHC-F 1:100-1:500 IHC-P 1:100-1:500
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20 C for long-term storage. Store at 2-8 C for up to one month.

TARGET INFORMATION


Gene ID	81876
Gene Symbol	RAB1B
UniProt ID	RAB1B_HUMAN
Immunogen Region	110aa C-Term
Immunogen Sequence	QEIDRYASENVNKLVLGNKS EIDRYASENVNKLVLGNKSD IDRYASENVNKLVLGNKSDL DRYASENVNKLVLGNKSDLT RYASENVNKLVLGNKSDLTT
Specificity	Detects endogenous levels of Rab1b protein by Western blot in whole cell lysates and transfected cells with GFP-Rab1b cds.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Anti RAB1B Goat Polyclonal antibody immunohistochemistry human STJ140139




Anti RAB1B Goat Polyclonal antibody immunohistochemistry human STJ140139



Anti RAB1B Goat Polyclonal antibody immunohistochemistry mouse STJ140139



Anti RAB1B Goat Polyclonal antibody immunohistochemistry human STJ140139



Anti RAB1B Goat Polyclonal antibody western blot hela STJ140139