

## Anti-PTK2 antibody (690aa C-Term) (STJ140154)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Canine/Monkey

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20 C for long-term storage. Store at 2-8 C for up to one month.

### TARGET INFORMATION

Gene ID	<a href="#">5747</a>
Gene Symbol	<a href="#">PTK2</a>
UniProt ID	<a href="#">FAK1_HUMAN</a>
Immunogen Region	690aa C-Term
Immunogen Sequence	MRMESRRQATVSWDSGGGSDE RMESRRQATVSWDSGGGSDEA MESRRQATVSWDSGGGSDEAP ESRRQATVSWDSGGGSDEAPP SRRQATVSWDSGGGSDEAPPK RRQATVSWDSGGGSDEAPPKP RQATVSWDSGGGSDEAPPKPS QATVSWDSGGGSDEAPPKPSR ATVSWDSGGGSDEAPPKPSRP TVSWDSGGGSDEAPPKPSRPG VSWDSGGGSDEAPPKPSRPGY SWDSGGGSDEAPPKPSRPGY
Specificity	Detects endogenous levels of PTK2 by Western blot.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti PTK2 Goat Polyclonal antibody western blot  
STJ140154