

Anti-TUBB antibody (345aa C-Term) (STJ140156)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Canine/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	203068
Gene Symbol	TUBB
UniProt ID	TBB5_HUMAN
Immunogen Region	345aa C-Term
Immunogen Sequence	MKEVDEQMLNVQKNSSYFV KEVDEQMLNVQKNSSYFVE EVDEQMLNVQKNSSYFVEW VDEQMLNVQKNSSYFVEWI DEQMLNVQKNSSYFVEWIP EQMLNVQKNSSYFVEWIPN QMLNVQKNSSYFVEWIPNN
Specificity	Detects endogenous levels of all variants (var1 52 KDa, var2 49.5 KDa, var3 27.3 KDa, var4 44.8KDa, var5 41.7 KDa, var6 41.7 KDa) of TUBB in different cell lines by Western blot.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti TUBB Goat Polyclonal antibody western blot
mcf7 STJ140156