

Anti-PDIA6 antibody (220aa C-Term) (STJ140160)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Canine/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	10130
Gene Symbol	PDIA6
UniProt ID	PDIA6_HUMAN
Immunogen Region	220aa C-Term
Immunogen Sequence	MDATVNQVLASRYGIRGFPT DATVNQVLASRYGIRGFPTI ATVNQVLASRYGIRGFPTIK TVNQVLASRYGIRGFPTIKI VNQVLASRYGIRGFPTIKIF NQVLASRYGIRGFPTIKIFQ QVLASRYGIRGFPTIKIFQK VLASRYGIRGFPTIKIFQKG LASRYGIRGFPTIKIFQKGE ASRYGIRGFPTIKIFQKGES SRYGIRGFPTIKIFQKGESP
Specificity	Detects endogenous levels of PDIA6 by Western blot.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti PDIA6 Goat Polyclonal antibody western blot
 hela STJ140160