

Anti-SNCB antibody (80 aa C-Term) (STJ140169)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Canine/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:5000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	6620
Gene Symbol	SNCB
UniProt ID	SYUB_HUMAN
Immunogen Region	80 aa C-Term
Immunogen Sequence	SRETGLVKREEFPTDLKPEE RETGLVKREEFPTDLKPEEV ETGLVKREEFPTDLKPEEVA
Specificity	Detects endogenous levels of SNCB by Western blot. This Ab is specific for SNCB and does not recognize SNCA and SNCG proteins.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Anti SNCB Goat Polyclonal antibody western blot mouse brain STJ140169