

Anti-BCL2 antibody (143 aa N-Term) (STJ140179)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Canine/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	596
Gene Symbol	BCL2
UniProt ID	BCL2_HUMAN
Immunogen Region	143 aa N-Term
Immunogen Sequence	MAHAGRTGYDNREIVMKYIH AHAGRTGYDNREIVMKYIHY HAGRTGYDNREIVMKYIHYK AGRTGYDNREIVMKYIHYKL GRTGYDNREIVMKYIHYKLS RTGYDNREIVMKYIHYKLSQ TGYDNREIVMKYIHYKLSQR GYDNREIVMKYIHYKLSQRG
Specificity	Detects a band of approximately 25 kDa by Western blot in the following cell lysate: LS174T, Jurkat and HUH.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti BCL2 Goat Polyclonal antibody western blot
STJ140179