

Anti-SARS-CoV-2 NSP2 antibody (480-stop) (STJ140183)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Region	480-stop
Immunogen Sequence	MGGQIVTCAKEIKESVQTFKLV GGQIVTCAKEIKESVQTFKLV GQIVTCAKEIKESVQTFKLV QIVTCAKEIKESVQTFKLV IVTCAKEIKESVQTFKLVN VTCAKEIKESVQTFKLVNK TCAKEIKESVQTFKLVNKF CAKEIKESVQTFKLVNKF
Specificity	In lysates of transfected cells with the plasmid containing the sequence used, detects the fusion protein by Western blot.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti SARS-CoV-2 NSP2 Goat Polyclonal antibody western blot STJ140183